
Elite2

Instructions for Use

EL24L1-EL30R9
EL24L1D-EL30R9D

EN	Instructions for Use	2
FR	Instructions D'Utilisation	14
DE	Gebrauchsanweisung	26
IT	Istruzioni per L'Uso	39
ES	Instrucciones de Uso	51

Blatchford:

Application

These instructions for use are for the practitioner.

The term *device* is used throughout these instructions for use to refer to Elite2.

This device is to be used exclusively as part of a lower limb prosthesis.

A high-energy-return foot. The independent heel and toe springs provide axial deflection. The split toe provides good ground compliance.

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Appropriate footwear must be worn to avoid the risk of slipping in wet environments.

Contra-indications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Intended for a single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

Spring set selection

		User weight										kg (lb)
Impact	Activity	44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-130 (256-285)	131-147 (286-325)	148-166 (326-365)	
Low	3	1	1	2	3	4	5	6	7	8	9	Foot spring set
Mod	4	1	2	3	4	5	6	7	8	9		
High	4	2	3	4	5	6	7	8	9			
Low	Daily walking and occasional sports such as golf and hiking											Important: For higher impact users, do not exceed the weight limit for individual springs.
Moderate	Aggressive walking, frequent or daily sports such as jogging											
High	Daily activities such as distance running, climbing, lifting and carrying heavy objects for vocational purposes											

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.


Note:

If in doubt choosing between two categories, choose the higher rate spring set.

Foot Spring set recommendations shown are for trans-tibial users.

For trans-femoral users we suggest selecting a spring set one category lower, refer to fitting advice Section 7 to ensure satisfactory function and range of movement


When fitted with springs, cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

2 Construction

Principal parts:

- Carrier Assembly (aluminum/St. Stl./titanium)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (titanium)
- Glide Sock (UHM PE)
- Foot Shell (PU)

3 Function

This device comprises an e-carbon toe and independent heel spring.

Heel and toe springs are attached to the carrier using titanium screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

4 Maintenance

Maintenance must be carried out by competent personnel.

We recommend the following annual maintenance:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Check heel and toe springs for signs of delamination or wear and replace if necessary. Some minor surface damage may occur after a period of use, this does not affect the function or strength of the foot.

The user must be handed the user information card supplied, and be advised of the following:

- Any changes in performance of this device must be reported to the practitioner e.g. reduced energy return or unusual noises.
- The practitioner must also be informed of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

5 Limitations on Use

Intended life

A local risk assessment should be carried out based upon activity and usage.

Lifting loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Foot units must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, it should be inverted and dried before further use.


Suitable for submersion

Exclusively for use between -15 °C and 50 °C.

We recommend using Blatchford products with this device.

6 Bench Alignment

Static alignment

Setup length

With flexion, adduction and abduction properly accommodated, set the limb length 5 mm longer than the sound side to allow for compression and deflection of the foot springs during gait. This should be re-assessed once the dynamic trial commences and the length adjusted accordingly.

Build line

This should fall 10 mm anterior to the centre of the pyramid (with heel height properly accommodated). The socket should be positioned accordingly.

Dynamic alignment

Coronal plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

Trans-femoral alignment

Align trans-femoral components according to fitting instructions supplied with the knee, keeping the build line relative to the foot as shown.

* Allow for users own footwear

7 Fitting Advice

Spring sets are supplied as matched pairs i.e. the heel and toe spring are designed to work together to give smooth progression for most users.

Heel wedge

A heel wedge is supplied with the foot. Fitting the wedge will have the effect of stiffening the heel spring. These can be taped in place for trial. For permanent fitting, wedges should be adhered in place by application of Loctite 424 (926104) between the lower contacting surface of the heel and the wedge.

Heel stiffness

Progression throughout the stance phase should be smooth; heel function is key to this process:


- Too soft a heel or load line excessively posterior will result in sinking at heel strike and difficulty in getting over the toe.
- Too hard a heel or load line excessively anterior will result in a rapid progression through mid-stance or jarring at heel strike.

	Symptoms	Remedy
Heel too soft	• Sinking at heel strike - Difficulty climbing over the toe (toe feels too hard)	1. Add heel wedge 2. Move socket anteriorly in relation to the foot (excess movement may result in drop off) 3. If 1 and 2 fail, fit a stiffer spring set
Heel too hard	• Rapid transition from heel strike through stance phase - Difficulty in controlling heel action, foot jars into mid-stance • Foot feels too rigid	1. Remove heel wedge (if fitted) 2. Move socket posteriorly in relation to foot 3. If 1 and 2 fail, fit softer spring set

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

8 Dis-assembly Instructions

1

Carefully remove, with a knife, any foam cosmesis that may be adhered to the foot shell.

2

Remove securing screws from the carrier and pull the toe spring towards the rear of the foot.

3

Remove the toe spring to leave the carrier/heel spring assembly alone inside the foot shell.

4

Rotate the carrier/heel spring assembly toward the rear of the foot to dis-engage the spring from it's location in the shell.


9 Assembly Instructions

1

Assemble heel spring onto carrier, use Loctite 243 (926012) and torque to 15 Nm. Use special Allen wrench 940080, Torque spanner adaptor 940081 or 13A/F spanner 940273.

2

Slide toe spring into position on carrier (do not assemble bolts). Fold the glide sock neatly around the end of the toe.

3

Unwrap glide sock over heel spring as shown.

4


If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

5

Push toe spring forward as far as possible.


6

Slide carrier/heel spring assembly forward into slot in top of foot shell.


9 Assembly Instructions *(continued)*

7


Slide a metal shoe horn between heel spring and back of foot shell opening and lever the spring down into the foot shell.

8


Push toe spring towards carrier. Ensure glide sock doesn't get trapped between the spring and carrier. Push spring into location in shell.

9

Press heel spring into location in shell as shown.

10


11

10 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	395 g (14 oz)
Activity Level:	3–4
Maximum User Weight:	166 kg (365 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Adjustment:	±7° Angular
Build Height: (See diagram below)	130 mm
Heel Height:	10 mm

Fitting length

11 Ordering Information

Order example

EL

25L

3

e.g. EL25L3

Size

Spring Set

Side

Category

Available from size 24 to size 30:
EL24L1 to EL30R9
EL24L1D to EL30R9D
(add 'D' for a dark tone foot shell)

Spring kit

Spring kit	Foot size		
	24–26	27–28	29–30
Set 1	539710	539719	Special order
Set 2	539711	539720	Special order
Set 3	539712	539721	539730
Set 4	539713	539722	539731
Set 5	539714	539723	539732
Set 6	539715	539724	539733
Set 7	539716	539725	539734
Set 8	539717	539726	539735
Set 9	539718	539727	539736

Foot shell

For dark add 'D'

Left	
24L	539005
25L	539007
26L	539009
27L	539011
28L	539013
29L	539015
30L	539017

Right	
24R	539006
25R	539008
26R	539010
27R	539012
28R	539014
29R	539016
30R	539018

Glide Sock

One size fits all Part No. 532811

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions. See Blatchford website for the current full warranty statement.

Trademark Acknowledgements

Elite Foot and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

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Blatchford SAS
Endolite India Ltd