
# EliteVT

## Instructions for Use

EVT24L11–EVT30R94  
EVT24L11D–EVT30R94D

|    |                            |    |
|----|----------------------------|----|
| EN | Instructions for Use       | 2  |
| FR | Instructions D'Utilisation | 14 |
| DE | Gebrauchsanweisung         | 26 |
| IT | Istruzioni per L'Uso       | 39 |
| ES | Instrucciones de Uso       | 51 |

**Blatchford:**

# 1 Description and Purpose

EN

## **Application**

These instructions are for use by the practitioner.

The term *device* is used throughout these instructions for use to refer to EliteVT.

This device is to be used exclusively as part of a lower limb prosthesis.

A high-energy-return foot. Independent heel and toe spring provides axial deflection. The split toe provides good ground compliance. This device also includes a spring element capable of absorbing axial and rotational shock intended to reduce sheer stresses at the residuum/socket interface.

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

## **Contra-indications**

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Intended for a single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

Spring set selection

|        |          | User weight        |                    |                    |                    |                    |                     |                      |                      |                      |                      |                                                                                   |                                        | Foot spring set | kg<br>(lbs) |
|--------|----------|--------------------|--------------------|--------------------|--------------------|--------------------|---------------------|----------------------|----------------------|----------------------|----------------------|-----------------------------------------------------------------------------------|----------------------------------------|-----------------|-------------|
| Impact | Activity | 44-52<br>(100-115) | 53-59<br>(116-130) | 60-68<br>(131-150) | 69-77<br>(151-170) | 78-88<br>(171-195) | 89-100<br>(196-220) | 101-116<br>(221-255) | 117-130<br>(256-285) | 131-147<br>(286-325) | 148-166<br>(326-365) |                                                                                   |                                        |                 |             |
| Low    | 3        | 1                  | 1                  | 2                  | 3                  | 4                  | 5                   | 6                    | 7                    | 8                    | 9                    |  | ◀ Axial shock spring<br>rated as shown |                 |             |
| Mod    | 4        | 1                  | 2                  | 3                  | 4                  | 5                  | 6                   | 7                    | 8                    | 9                    |                      |                                                                                   |                                        |                 |             |
| High   | 4        | 2                  | 3                  | 4                  | 5                  | 6                  | 7                   | 8                    | 9                    |                      |                      |                                                                                   |                                        |                 |             |

Low                    Daily walking and occasional sports such as golf and hiking

Moderate           Aggressive walking, frequent or daily sports such as jogging

High                  Daily activities such as distance running, climbing, lifting and carrying heavy objects for vocational purposes

**Important: For higher impact users, do not exceed the weight limit for individual springs.**

**Note:**

If in doubt choosing between two categories, choose the higher rate spring set.

Foot spring set recommendations shown are for trans-tibial users. For trans-femoral users we suggest selecting a foot spring set one category lower.


Once fitted with springs, cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

## 2 Construction

### Principal parts:

- Carrier Assembly (aluminum/ titanium/ St. Stl.)
- Heel & Toe Springs (e-carbon)
- Spring Attachment Screws (titanium)
- Glide Sock (UHM PE)
- Foot Shell (PU)

---

## 3 Function

This device comprises an e-carbon toe and independent heel spring. The carrier assembly contains a spring component which allows the pyramid tube to deflect and/or rotate through a limited range under an axial and/or torsional load. When the load is removed the spring returns to its original position. The effect of this movement is to attenuate any shock forces that would otherwise be transmitted to the residuum/socket interface both in axial and rotational directions. Heel and toe springs are attached to the carrier assembly using titanium screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

As part of routine maintenance to be carried out annually as a minimum, the following is recommended:

- Remove foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check heel and toe springs for signs of damage or wear.

Re-lubrication of the bearing surface may be carried out using 'Sapphire Endure' (928015) or equivalent PTFE loaded grease, if necessary, by disassembling and reassembling as shown in the Axial Shock spring replacement instructions.

The user should be advised:-

Any changes in performance of this device must be reported to the practitioner e.g. unusual noises, increased stiffness or restricted/excessive rotation.

The practitioner must also be informed of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

---

## 5 Limitations on Use

### Intended life

A local risk assessment should be carried out based upon activity and usage.

### Lifting loads

Users weight and activity is governed by the stated limits.

Load carrying by the users should be based on a local risk assessment.

### Environment

Avoid exposing this device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

We recommend using Blatchford products with this device.

Suitable for outdoor use

## 6 Bench Alignment

## Static alignment

Setup length

With flexion, adduction and abduction properly accommodated, set the limb length 10 mm longer than the sound side to allow for compression and deflection of the foot & axial shock springs during gait. This should be re-assessed before the dynamic trial commences and the length adjusted accordingly.

## Build line

This should fall 10 mm anterior to the centre of the pyramid (with heel height properly accommodated). The socket should be positioned accordingly.

## Dynamic alignment

## Coronal plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

## Sagittal plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

## Trans-femoral alignment

Align trans-femoral components according to fitting instructions supplied with the knee, keeping the build line relative to this device, as shown.

- \* Allow for users own footwear

# 7 Fitting Advice

Springs are supplied as matched sets i.e. the heel, toe and axial springs are designed to work together to give smooth progression for most users.

## Heel Wedge

A heel wedge is supplied with the foot. Fitting the wedge will have the effect of stiffening the heel spring. These can be taped in place for trial. For permanent fitting, wedges should be adhered in place by application of Loctite 424 (926104) between the lower contacting surface of the heel and the wedge.

## Axial Stiffness

The axial spring controls both compressive and torsional resilience. Typical vertical movement should be between 5–10 mm for normal walking. Alternative stiffness springs are available if required.

| Condition     | Symptoms                                                                                                                                           | Remedy                                                                                                                                                          |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Heel too soft | Sinking at heel strike<br>Difficulty climbing over the toe (toe feels too hard)                                                                    | 1. Add heel wedge<br>2. Move socket anteriorly in relation to the foot (excess movement may result in drop-off)<br>3. If 1 and 2 fail, fit a stiffer spring set |
| Heel too hard | Rapid transition from heel strike through stance phase<br>Difficulty in controlling heel action, foot jars into mid-stance<br>Foot feels too rigid | 1. Remove heel wedge (if fitted)<br>2. Move socket posteriorly in relation to foot<br>3. If 1 and 2 fail, fit softer spring set                                 |

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

## 8 Foot Shell Removal & Foot Spring Replacement

1


Remove toe spring screws (13A/F spanner 940273).

2


Withdraw toe spring from shell.

3

Rotate the carrier/heel spring assembly toward the rear of the foot to dis-engage the spring from it's location in the shell.

4

Undo and remove heel spring screws. Assemble new heel spring onto carrier, use Loctite 243 (926012) and torque to 15 Nm. Use special Allen wrench 940080, Torque spanner adaptor 940081 or 13A/F spanner 940273

5

Loosely assemble toe spring and glide sock, mark carrier with spring category as appropriate.

6

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.


## 8 Foot Spring Replacement (continued)

7


Insert assembly into foot shell, push toe spring as far forward as possible while sliding carrier and heel assembly forward in the foot shell.

8

Slide a metal shoe horn between heel spring and back of foot shell opening and lever the spring down into the foot shell.


9

Toe Spring location in Foot Shell

Push toe spring towards carrier. Ensure glide sock doesn't get trapped between the spring and carrier. Push spring into location in shell.

10

Heel Spring location slot

Press heel spring into location in shell as shown.

11

Attach the toe spring to the carrier using Loctite 243 (926012) on the bolts.

1. For spring rates 1 to 7 use 4 A/F Allen wrench and torque to 15 Nm. Do not use external hex, this is reserved for loosening the bolt, if required.
2. For spring rates 8 and 9 use 13 A/F spanner and torque to 25 Nm.





12

Bond foam cosmesis to top surface of the foot shell as shown using Thixofix adhesive (926204) or equivalent and shape to suit.

# 9 Axial Shock Spring Replacement

Remove from foot shell as shown previously


Re-lubricate bearing surface with 'Sapphire Endure' (928015) or equivalent PTFE loaded grease if necessary.

Re-assembly is the reverse of above, note the position and location of anti-rotation pins in the axial spring and the orientation of the slot in the pylon (posterior). Apply Loctite and apply torque to screws as shown on page 4.

# 10 Technical Data

|                                          |                                     |
|------------------------------------------|-------------------------------------|
| Operating and Storage Temperature Range: | -15 °C to 50 °C<br>(5 °F to 122 °F) |
| Component Weight (size 26):              | 1050 g ( 2lb 5 oz)                  |
| Activity Level:                          | 3–4                                 |
| Maximum User Weight:                     | 166 kg (275 lb)                     |
| Range of Adjustment                      | ±7° angular                         |
| Proximal Attachment:                     | Male pyramid (Blatchford)           |
| Build Height:                            | 170 mm                              |
| Heel Height                              | 10 mm                               |

## Fitting length

\* sizes  
22–28 = 70 mm  
29–30 = 80 mm

# 11 Ordering Information

## Order example

EVT

25L

7

3

Size  
Side

Spring  
Set

Shock  
spring

e.g. EVT25L73

Available from size 24 to size 30:  
EVT24L11 to EVT30R94  
EVT24L11D to EVT30R94D  
*(add 'D' for a dark tone foot shell)*

## Spring Kit

| Spring Kit | Foot Size |        |               |
|------------|-----------|--------|---------------|
|            | 24–26     | 27–28  | 29–30         |
| Set 1      | 539710    | 539719 | Special order |
| Set 2      | 539711    | 539720 | Special order |
| Set 3      | 539712    | 539721 | 539730        |
| Set 4      | 539713    | 539722 | 539731        |
| Set 5      | 539714    | 539723 | 539732        |
| Set 6      | 539715    | 539724 | 539733        |
| Set 7      | 539716    | 539725 | 539734        |
| Set 8      | 539717    | 539726 | 539735        |
| Set 9      | 539718    | 539727 | 539736        |

## Foot Shell

For dark add 'D'

| Left |        |
|------|--------|
| 24L  | 539005 |
| 25L  | 539007 |
| 26L  | 539009 |
| 27L  | 539011 |
| 28L  | 539013 |
| 29L  | 539015 |
| 30L  | 539017 |

| Right |        |
|-------|--------|
| 24R   | 539006 |
| 25R   | 539008 |
| 26R   | 539010 |
| 27R   | 539012 |
| 28R   | 539014 |
| 29R   | 539016 |
| 30R   | 539018 |

## Axial Shock Springs

|        |        |
|--------|--------|
| Rate 1 | 539058 |
| Rate 2 | 539059 |
| Rate 3 | 539060 |
| Rate 4 | 539061 |

## Glide Sock

|                   |        |
|-------------------|--------|
| One size fits all | 532811 |
|-------------------|--------|

**Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

**CE conformity**

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class I Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

**Warranty**

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

**Trademark Acknowledgements**

Elite Foot and Blatchford are registered trademarks of Blatchford Products Limited.

**Manufacturer's Registered Address**

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

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Blatchford SAS
Endolite India Ltd