Javelin

Instructions for Use

JAV22L1S–JAV30R9S
JAV22L1SD–JAV30R9SD

EN	Instructions for Use	2
FR	Instructions D'Utilisation	14
DE	Gebrauchsanweisung	26
IT	Istruzioni per L'Uso	39
ES	Instrucciones de Uso	51

1 Description and Purpose

Application

These instructions are for the practitioner.

The term *device* is used throughout these instructions for use to refer to Javelin.

This device is to be used exclusively as part of a lower limb prosthesis.

A lightweight, high-energy-return foot, ideal for medium impact use. The flexible extended pylon reduces shock load transmission making it suitable for work and leisure activities. Independent heel and toe springs provide some axial deflection. The split toe provides good ground compliance.

This device is recommended for users that have the potential to achieve Activity Level 3. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in the Activity Levels 2 and 4* who would benefit from the smoother transition from heel to toe offered by this device, but this decision should be made with sound and thorough justification.

*Maximum user weight 100 kg and always use one higher spring rate category than shown in the Spring set selection table.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

Contra-indications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Intended for a single user.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

Spring set selection

		User weight										
Activity		44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166	kg
		(100-115)	(116-130)	(131-150)	(151-170)	(171-195)	(196-220)	(221-255)	(256-285)	(286-325)	(326-365)	(lb)
3		1	1	2	3	4	5	6	7	8	9	Foot Spring set

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Note:

If in doubt choosing between two categories, choose the higher rate spring set.


Foot Spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a spring set one category lower, refer to fitting advice Section 9 to ensure satisfactory function and range of movement

2 Construction

Principal parts:

- Carrier (aluminum)
- Heel & Toe Springs (e-carbon)
- Spring attachment screws (titanium)
- Glide Sock (UHM PE)
- Foot shell (PU)

3 Function

This device comprises an e-carbon toe spring pylon and independent heel spring. The pylon can be cut to length to suit individual users.

Heel and toe springs are attached to the carrier using titanium screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

4 Maintenance

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use, this does not affect the function or strength of the foot.

The user must be handed the user information card supplied, and be advised of the following:

- Any changes in performance of this device must be reported to the practitioner e.g. reduced energy return or unusual noises.
- The practitioner must also be informed of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

5 Limitations on use:

Intended life

A local risk assessment should be carried out based upon activity and usage.

Lifting loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Foot products must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, it should be inverted and dried before further use.


Suitable for submersion

Exclusively for use between -15 °C and 50 °C (5 °F and 122 °F).

We recommend using only Blatchford products with this device.

6 Bench Alignment

Static alignment

Setup length

Note: Shin pylon should be vertical when used with a 10 mm heel height shoe. Lower heel height shoes may require the use of a packing piece below the heel.

Cut shin to approximate length, allowing 10–15 mm for final trimming. Note trial length will typically include 5 mm to allow for axial compression of heel and toe.

Build line

This should fall 1/3 of the foot length from the heel.

Dynamic alignment

Coronal plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

Transfemoral alignment

Align transfemoral components according to fitting instructions supplied with the knee, keeping the build line relative to the pylon as shown.

* Allow for users own footwear

7 Cutting Pylon to Length

Note: The following instructions assume the use of the male, sliding, rotating pyramid. If an alternative adaptor is required, this should be taken into account throughout the build process.

Use appropriate health and safety equipment at all times including extraction facilities

1

 Cut longer than required

 5

Cut the pylon 10–15 mm longer than required and assemble fitting clamp (941255).

2

 4

15 Nm

Fit socket to fitting clamp top housing and tighten pyramid screws (15 Nm). Refer to alignment page for initial set up.

3

 5

15 Nm


Tighten fitting clamp bolts to 15 Nm (check each screw at least twice).

4

Stand user, carry out dynamic alignment and walking trials. Adjust height to achieve definitive length with weight evenly distributed.

5

Remove socket, place fitting clamp in the vice and cut blade to length. Fit shoulder screws to prevent damage to components


6

Drill two 8.1 dia. holes through drill bushes.

7 Cutting Pylon to Length *(continued)*

7

Deburr end of pylon.

8

Slide pyramid into top housing (or fit preferred alignment option).

9

4
25 Nm

Loctite (243) & tighten (25 Nm) top housing bolts. Attach and align socket.

8 Foot Shell Removal & Pylon/Spring Replacement

Be aware of finger trap hazard at all times.

1

4
Remove top housing from pylon.

2

4
Remove lower carrier bolt.

8 Foot Shell Removal & Pylon/Spring Replacement *(continued)*

3

Remove pylon from carrier/foot shell.

4

Rotate heel/carrier assembly from foot shell/glide sock.

5

Remove heel spring from carrier.

6

Assemble replacement heel spring, applying Loctite 243 to screw threads and torque to 15 Nm.

7




Attach pylon to carrier. Loctite (243) & tighten bolt (25 Nm).

8

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

8 Foot Shell Removal & Pylon/Spring Replacement *(continued)*

9 Fitting Advice

Springs are supplied as matched sets i.e. the pylon/toe and heel springs are designed to work together to give smooth progression for most users. If this is not the case, please refer to the advice below.


	Symptoms	Remedy
Heel too soft	<ul style="list-style-type: none">• Sinking at heel strike• Difficulty climbing over the toe (toe feels too hard)	Move socket anteriorly in relation to the foot (excess movement may result in drop-off)
Heel too hard	<ul style="list-style-type: none">• Rapid transition from heel strike through stance phase• Difficulty in controlling heel action, foot jars into mid-stance• Foot feels too rigid	Move socket posteriorly in relation to foot
Toe/Pylon too soft	<ul style="list-style-type: none">• Rapid progression through mid-stance.• <i>Drop-off</i> at higher activity levels	1 Move socket posteriorly in relation to foot 2 Plantarflex foot slightly Note: some re-alignment may be required

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

10 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (Size 26N):	670 g (1 lb 8 oz)
Activity Level:	3
Maximum User Weight:	166 kg (365 lb)
Proximal Alignment Attachment:	Available with separately supplied alignment devices
Build Height: (See diagram below)	200–365 mm
Heel Height	10 mm

Fitting Length

Size	A
22–26	65 mm
27–28	70 mm
29–30	75 mm

11 Ordering Information

Order example

JAV	25	L	N	3	S
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

JAV22L1S to JAV30R9S

JAV22L1SD to JAV30R9SD

(add 'D' for a dark tone foot shell)

*Sizes 25-27 only. For all other sizes, omit the Width field.

e.g. JAV25LN3S, JAV22R4S, JAV27RW4SD

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22-24	25-26	27-28	29-30
Set 1	539301S	539510S	539519S	Special Order
Set 2	539302S	539511S	539520S	Special Order
Set 3	539303S	539512S	539521S	539530S
Set 4	539304S	539513S	539522S	539531S
Set 5	539305S	539514S	539523S	539532S
Set 6	539306S	539515S	539524S	539533S
Set 7	539307S	539516S	539525S	539534S
Set 8	539308S	539517S	539526S	539535S
Set 9	-	539518S	539527S	539536S

Foot Shell (for dark add 'D')		
Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	-	539050S
28R	-	539051S
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

Item	Part. No.
Glide Sock (Sizes 22-26)	531011
Glide Sock (Sizes 27-30)	532811
Fitting Clamp	941255

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Manufacturer's Registered Address

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Blatchford SAS
Endolite India Ltd