# Elan

## Clinician's Manual

ELAN22L1S—ELAN30R8S  
ELAN22L1SD—ELAN30R8SD

EN	Clinician's Manual	1
FR	Manuel de L'Orthoprothésiste	33
DE	Gebrauchsanweisung für Techniker	65
IT	Manuale per il Tecnico Ortopedico	97
ES	Manual del Técnico	129
NO	Klinikerhåndbok	161
RU	Руководство протезиста	193
TR	Klinisyen Kılavuzu	225

# 1 Description and purpose

EN

The Elan is a biomimetic foot and ankle system with microprocessor controlled dorsiflexion and plantar flexion. It is designed to adapt dynamically (via motorized valves) to provide assistance when walking on a variety of surfaces, when changing gait speeds and increased levels of safety and security when walking and standing.

## Application

- These instructions are for use by the practitioner
- The Elan is to be used exclusively as part of a lower limb prosthesis
- Intended for single user

Ensure that the user has understood all instructions for use, drawing particular attention to the sections regarding maintenance and use of the battery charger connection.

## Features

*(See also Section 3 Function)*

- Dynamic damping resistance aims to improve gait characteristics enabling safer walking
- Microprocessor control of resistance in real time
- The Elan Software adapts to individual user gait characteristics
- Standing Support mode
- Aims to reduce stress at the socket interface
- PC interface setup with Bluetooth® link and direct communications
- Power is supplied via an integral rechargeable battery pack

This device provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. It is intended to improve postural sway, balance stability and inter-limb loading symmetry while easing high pressures at the socket interface.

## Activity level

This device is recommended for users that have the potential to achieve Activity Level 3 who may benefit from enhanced stability and an increase in confidence on uneven surfaces.

There are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users of Activity Levels 2 and 4\* who would benefit from the enhanced stability offered by the Elan.

Users at Activity Levels 2 and 4\* will require softer or stiffer springs as appropriate for the individual user rather than as shown in the spring selection guide.

\*(maximum user weight 100 kg (220 lbs) and always use one higher spring rate category than shown in the Spring Set Selection table)

## Contra-indications

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Activity	User Weight*								kg (lbs)	Foot Spring set
	44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-125 (256-275)		
<b>3</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>		

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Note:

In doubt choosing between two categories, choose the higher rate spring set.

Foot spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a foot spring set one category lower, refer to fitting advice Section 9.6 to ensure satisfactory function and range of movement.

### Order Example:

<b>ELAN</b>	<b>25</b>	<b>L</b>	<b>N</b>	<b>3</b>	<b>S</b>
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:  
ELAN22L1S to ELAN30R8S  
ELAN22L1SD to ELAN30R8SD

\*Sizes 25-27 only. For all other sizes, omit the Width field.

(add 'D' for a dark tone foot shell)

e.g. ELAN25LN3S, ELAN22R4S, ELAN27RW4SD

### Checking the package contents:

1. Elan
2. Battery Charger and AC adapters
3. Charger Cable adapter
4. Clinician's Manual
5. User Guide
6. Bluetooth Dongle (Micro-D)
7. x2 Micro-D Connector Cover
8. External USB Bluetooth® Adapter

Part No.  
409087E

### Also required for Elan set-up:

(supplied separately Part No. **406431E**)

1. USB stick and interface software
2. Comms Cable (USB)
3. Comms Cable Adapter
4. Interface cable (Micro-D)
5. External USB Bluetooth® Adapter

\* refer to Section 11 Technical Data

# Safety Information


-  **The caution symbol highlights safety information which must be followed carefully.**
-  1. Any changes in performance of the ankle motion e.g. restricted movement, non-smooth motion should be immediately reported to your service provider.
  -  2. Always use a hand rail when descending stairs and at any other time if available.
  -  3. Any excessive changes in heel height after programming will adversely affect ankle function and should be immediately reported to your service provider for reprogramming and calibration.
  -  4. After continuous use the ankle casing may become hot to the touch.
  -  5. Avoid strong magnetic fields, sources of electrical interference, atmospheres containing liquids and/or powders.
  -  6. Avoid exposure to extreme heat and/or cold.
  -  7. The electronic components are sealed for water resistance. However the Elan is not intended for use when immersed in water or as a shower prosthesis. If the ankle comes into contact with water wipe dry immediately.
  -  8. Ensure no water enters the interface/charging cable connector.
  -  9. The Elan is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.
  -  10. Contact your service provider for advice if the Elan is to be used with a pacemaker or any other electronic medical device.
  -  11. Assembly, programming and repair of the Elan must only be carried out by suitably qualified practitioners, certified by Blatchford.
  -  12. The user must not adjust or tamper with the setup of the Elan.
  -  13. The user should be advised to contact their clinician if their condition changes.
  -  14. The Elan detects standing and walking on inclines, declines and walking at various speeds on normal ground. Non-identified movements such as treadmills and travalators may cause unexpected ankle behaviour.
  -  15. Never plug the charger into the Elan while still wearing the limb.
  -  16. Only charge the batteries in the Elan using the charger supplied with the product. Do not use this charger for charging any other device. Ensure the charging port connector cover is correctly refitted after charging, replace if damaged.
  -  17. Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.
  -  18. Ensure charging cable and Bluetooth dongle are correctly aligned prior to connection and that adequate access to the charging port is provided (See Section 6).
  -  19. To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.
  -  20. Be aware of finger trap hazard at all times

## 2 Construction

### Principal parts:

- Hydraulic Body Assembly including pyramid (aluminum/St. Stl./titanium incorporates batteries (NiMH), electric motors and PCB
- Carrier Assembly (aluminum/St. Stl.)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (titanium/St. Stl.)
- Glide Sock (UHM PE)
- Foot Shell (PU)

### Component Identification

### 3 Function

With microprocessor control the Elan ankle adapts dynamically (via motorised valves) the ankle plantar flexion and dorsiflexion resistance settings to match the user's walking mode and terrain.

Activity	Elan Mode	Hydraulic Resistance
Walking (level ground)	Additional ground clearance in swing phase.	Slightly dorsiflexed Basic resistance settings
Walking fast	<i>Fast Walk</i> <i>[Assist]</i> Added momentum in the 'step to step' transition of the gait cycle ( <i>transtibial only</i> ).	Increased plantar flexion resistance, decreased dorsiflexion resistance
Walking uphill (moderate/steep incline)	<i>Ramp Assist</i> Enables easier walking up inclines.	Progressively changes; increasing plantar flexion and decreasing dorsiflexion resistance
Walking downhill	<i>Ramp Brake</i> Helps forward momentum and provides greater stability and security	Progressively changes; decreasing plantar flexion resistance and increasing dorsiflexion resistance
Standing	Basic Resistance settings or <i>Standing Support*</i> assistance	User settings from set-up/ high plantar flexion and dorsiflexion resistance
Low battery power	<i>Sleep</i>	Reverts to Basic Resistance settings from set-up (No dynamic control)

Note:

Standing Support mode provides trip and fall risk reduction and prevention features.

\*This feature is only available when enabled in the software setup.

To enable Standing Support mode see Section 14 *Options*> *Standing Support Mode*.

---

## 4 Maintenance

Maintenance must be carried out by competent personnel.

A service is required every 20 months as a condition of warranty.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use; this does not affect the function or strength of the foot.

(See Section 10 for the foot shell removal and spring replacement procedure)

The wearer should be advised:

Any changes in performance of this device must be reported to the practitioner

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers. Ensure water/liquids do not enter the interface cable connector.

---

## 5 Limitations on use

### Intended life:

A local risk assessment should be carried out based upon activity and usage.

### Lifting loads:

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment:

Avoid exposing the Elan to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C (5 °-122 °F).

It is recommended that only Blatchford products be used in conjunction with the Elan.

## 6 Battery Charging


**i** Please advise user of the charging information in this section. The batteries are not field replaceable and will be replaced as necessary when serviced.

### Read the user instructions supplied with the charger before using for the first time.

We recommend the Elan is charged daily; it should take 2 hours to charge fully, depending on residual charge. An exhausted or flat battery may take longer.

Lay the limb on a flat surface and use the Charger Cable Adapter to connect the charger into the socket in the top of the Elan body.

Use the indicator light on the charger to confirm that charging has initiated correctly:

LED	Charging state
Yellow	Initialisation state/no battery
Orange	Fast charge
Green/Yellow	Maintain/ Trickle charge
Green	Fully charged
Orange/Green	Error

The Elan will take several seconds to reset when the charger is disconnected and a single blue flash from the LED will confirm a correct startup.


- !** 1. Only charge the batteries in the Elan using the charger supplied with the product. Do not use this charger for charging any other device.
- !** 2. The battery charger can be used with an electrical outlet that is rated between 100V and 240V AC, 50/60Hz.
- !** 3. Never plug the charger into the Elan while still wearing the limb.
- !** 4. During charging no power is delivered to the ankle control system.
- !** 5. The front of the ankle may become warm during charging but should never be too hot to touch. There should be no abnormal smell. If this is the case immediately unplug from electrical supply and contact your service provider.
- !** 6. Do not refit the limb until the charger has been disconnected.
- !** 7. Always recharge the Elan before use if it has not been used for an extended period of time.

- ⚠ 8. Ensure the connector cover is refitted after charging
- ⚠ 9. When charging an empty battery the charger may revert to trickle charge as a precaution to prevent damaging the battery.

## Battery Charge/Status Indication

To check the battery charge status, tilt the foot down vertically and hold stationary for 4 seconds. The LED at the rear of the Elan will flash blue, 1-3 times (see table) and then repeat to indicate the charge status. Tilt the foot again to repeat the battery status check.

Note: this feature can be switched off, see Section 14 > Options.

LED	Battery Charge
3 flashes	Maximum
2 flashes	High
1 flash	Low

After charging the Elan the LED may indicate 2 or 3 flashes this will depend on the charging time, the age and condition of the battery.

## Elan Status Indication

### 1 Single Flash

When the charger is disconnected the LED will flash once to indicate correct start up.

### 2 Continuous Flash


If the LED flashes continuously after disconnection this may indicate an error has been detected. To reset the Elan, reconnect and then disconnect the charger the Elan should flash once to indicate correct start up. If the Elan will not reset contact your service provider.

### 3 Alternate Long-Short Flash

Indicates Battery shutdown; If this is seen shortly after an adequate charging period, contact your service provider a battery service may be required.

# 7 Bench Alignment


Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

Align transfemoral devices according to the instructions for use that are supplied with the knee.

### Tilt setting

Align limb to achieve range of motion shown.

$*$  Allow for users own footwear

## 8 Install the Endolite PC Software Interface

Note: when installing a newer version of the interface first uninstall any older versions.

To install the Endolite PC Software Interface on a computer running Microsoft Windows (see System requirements below) follow this procedure:

- 1 Put the PC software USB stick in a spare USB port.
- 2 In the folder *EN*

double click *setup.exe*

the *endolite Elan Setup Wizard* screen will display

click [*Next*]

and follow the on-screen message prompts to continue the installation.

- 3 When *Installation Complete* displays

click [*Close*]

to exit the install program.

### System requirements:

Bluetooth® 2.0 or later

Supported Operating Systems:

- Windows 8
- Windows 7
- SP1 Windows 7

Supported Architectures: x86, x64

Hardware Requirements:



Recommended Minimum: Intel Dual Core Duo

1 GHz or higher with 512 MB RAM or more

Minimum disk space: 100 MB

Prerequisites:

Windows Installer 3.1 or later  
Internet Explorer 5.01 or later

---

## 9 Using the PC Software Interface to set up the Elan

### Getting started:


The Elan is a biomimetic foot and ankle system with microprocessor control for plantar and dorsiflexion resistance. This device is designed to adapt dynamically to provide assistance when walking and standing on a variety of surfaces and changing gait speeds.

To set up the ankle resistance setting to an individual user preference the Elan control system is adjusted and calibrated with the PC interface via a Bluetooth® link.

The final walking calibration routine is used by the computer program in the Elan to adapt to gait characteristics after the alignment and resistance settings have been optimized. This calibration data is then used by the Elan software to recognize and adapt the Elan automatically to changing walking and daily living patterns.


### Set-up procedure:

- 1 Ensure the Elan batteries are charged (see Section 6).
- 2 Complete the Bench Alignment for the Elan (see Section 7).
- 3 Install the PC software interface (see Section 8) on a Microsoft Windows computer or remote device with a Bluetooth® capability.  
Ensure Bluetooth is enabled on the computer/remote device and disconnected to any other device.
- 4 Connect the Elan Bluetooth® Dongle to the Elan and ensure the blue LED is flashing ready to connect.  
Pair the Elan Dongle to your computer/remote device using the Bluetooth® manager program on your computer.  
If prompted during installation enter pair code > **1234**.  
Or you may be asked to confirm an ID code of the Elan Bluetooth dongle when connected > Yes  
After pairing disconnect the Dongle.  
If Bluetooth connection is not possible use a direct connection with the Comms Cable (USB). See procedure on the following page.
- 5 Carry out the Elan Alignment, Set-up and Calibration procedure using the PC interface controls with the user walking on a horizontal and level surface (See procedure on the following pages).

## Bluetooth Range


The Bluetooth range of operation may vary depending on the remote device used, we recommend a range of 5-10 meters distance. If the connection terminates unexpectedly and the LED indicator changes from GREEN (connected) to BLUE (standby), reduce the distance between the Elan and the remote device being used until the connection is restored and can be maintained.

## No Bluetooth® connection: using the Comms Cable (USB)

The Comms Cable can be used for setting up the Elan via a USB port with a laptop computer with no Bluetooth capability or where it is not possible to make a Bluetooth connection.

See Section 13 for the software driver installation procedure and instructions for using the Comms Cable.

## 9.1 Launch the PC Interface Program

1 To launch the PC Interface program:

Click the Endolite program icon in the list of programs *or* from your desktop.


2 When the *Login* screen appears

Click on: **B**

**i** Before setting up the Elan ensure that the battery is charged see section 6.

**A** Name and user Password are only required where access is restricted.

to continue to the *Start* screen.

## 9.2 Bluetooth® Connection to the Elan [Start]

1 Ensure the Elan is charged up and the Bluetooth® Dongle is connected to the Elan and flashing Blue\*.

2 Select the paired port for the Elan Dongle from the pull down menu. **A**

If this dongle has been used previously you can use [Autoconnect] **B**


When connected, the Dongle LED will change to green and the Bluetooth symbol on the status bar will change to active.

**i** To return to this screen click on:

**C** 

**i** Where it is not possible to make a Bluetooth connection use the Comms Cable via a USB port on the laptop.

See Section 13 for software driver installation procedure and instructions for using the Comms Cable.

### \*Note

If you have just paired the Bluetooth Dongle it may be necessary on some remote devices to disconnect and reconnect the Bluetooth Dongle to the Elan before making the Bluetooth connection from the software interface.

**D** [*Disconnect*]  
To disconnect the Bluetooth® Dongle from the PC/Laptop

**E** [*System Reset*]  
Resets the software in the Elan if the Bluetooth® connection is interrupted

## Status Bar

The current status of the Elan is shown by the indicators on the Status bar.

Mode Indicator	Description
Mode	Displays the current Elan mode
Angle	Shows foot angle dynamically
Tilt	Shows tilt position dynamically
PF and DF Valve	Displays the current dynamic resistance settings for plantar flexion and dorsiflexion
Battery 	shows the Elan is powered up
Power Save/ Sleep 	If displayed the Elan is in Power Save mode*
Standing Support* 	If displayed Standing Support Mode is active, high ankle resistance is selected.
Bluetooth 	Shows that the Elan is connected/disconnected via Bluetooth

### Power Save/ Sleep

After Calibration if the Elan detects no movement for approximately 10 seconds it activates power saving mode disabling the sensors. Note during the calibration and set-up the time interval is increased to 90 seconds.

### \*Standing Support Mode

This mode is only available if enabled, the default setting is OFF see Section 14 *Options*. When the Elan detects a stationary standing posture high ankle resistance is selected.

---

### To adjust or fine tune the setting on a previously calibrated Elan:

If the alignment has not been changed it is not necessary to follow the whole set up procedure again. Follow the **Walking Calibration** sequence (see section 9.5) and the Elan software will then add the new information to the previous data and average out the result giving improved settings for the user.

If alignment adjustments are required carry out a full factory reset and follow the full set up procedure.


## 9.3 Static Alignment

1 Complete the Bench Alignment set-up (see Section 7 )

2 Click on

to select the *Dynamic Alignment Vector* screen.

### 3 Static Alignment Procedure

With the user standing comfortably adjust the Elan A-P shift to achieve a neutral balance position\*.

See following pages.

#### Tilt **B**

Use the display to adjust the angular position of the Elan when standing.

Adjust using shift and/or tilt devices as required.


( When standing the *Vector* indicator position shown is recommended)

**Warning:**

Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

\* Ensure the Elan is uncalibrated prior to static alignment. See Section 9.5 Factory Reset.


#### Tilt setting position

# Biomimetic Alignment

The aim of alignment is to achieve a *balance point* while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved.

Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.


\* Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

# Biomimetic Adjustment

Note: Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.

Check the Bench Alignment and the heel height


Use shift for static alignment and standing. The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

## 9.4 Dynamic Alignment

Click on:

**A** Resistance Setup

to select the **Basic Resistance Settings** screen.

### Dynamic Alignment Procedure

1. Instruct the user to walk for approximately 3 minutes to become familiar with the resistance settings.
2. Adjust the Plantar flexion and Dorsiflexion hydraulic resistance settings to optimize the Elan to suit the user's preference for normal walking. [ Note the scale is from 4 to 7 ] **B**
3. Allow some time for practice and familiarization with the ankle resistance settings.
4. Verify that the hydraulics are moving when walking using the Vector Indicator. Reduce the resistance settings if minimal movement is shown. **C**


#### Power Save/ Sleep

If the Elan detects no movement for approximately 90 seconds it activates power saving mode, the Power Save symbol appears on the Status Bar and all the indicators will stop moving. Moving the Elan will automatically bring it back to normal operation.

## Dynamic Alignment:

### Resistance setting troubleshooting

## 9.5 Walking Calibration

Click on:

to select the *Walking Calibration* screen.

Note:

For earlier versions of the Elan the screen will be slightly different See Section 15.

 Calibration may be repeated at any time to fine tune the set-up and is recommended after any change in alignment and/or basic resistance settings.

### Walking Calibration procedure:


1 Move the slider to select the amputation level,

 Transfemoral *TF* or Transtibial *TT*

2 Start Calibration

To calibrate the Elan to the user's individual gait characteristics instruct the user to walk on a horizontal, level surface and at a comfortable pace.

After several steps, click on:

The step counter counts down the number of remaining steps required. 

(Note: 8-10 uninterrupted, regular steps on level ground are required for calibration).

 The Elan adapts to the individual user's gait characteristics from the Calibration stages. (See Section 3 for more information)

#### Factory Reset:

Should it be necessary to begin a completely new set-up for the Elan,

#### Caution:

Clicking on [**Factory Reset**] will **delete** all the Dynamic Alignment adjustment and any walking calibration data stored in the Elan.

 Click on 

If this option is chosen return to Section 9.4 to begin set-up.

## 9.6 Fitting advice

The correct alignment (A-P position), range of motion (distribution of plantar to dorsiflexion) and adjustment of the hydraulic resistance settings are critical in achieving a smooth roll-over and correct gradient adaptation.

The Elan foot is supplied assembled with heel and toe springs of the same category. If after following the instructions below you have problems with the function please contact the sales team in your area for advice.

Any of the following will have a negative effect on function and stability.

- Incorrect spring selection
- Incorrect A-P shift alignment
- Incorrect distribution of plantar and dorsiflexion range

Prior to setup or if further adjustment is required ensure the ankle is not calibrated. Ensure the ankle is not calibrated by carrying out a factory reset and that standing support mode is disabled. See Section 9.5 for details of factory reset and Section 14 to switch off standing support mode.


	Symptoms	Remedy
1	<ul style="list-style-type: none"><li>• Sinking at heel strike</li><li>• Difficulty in achieving a smooth progression to mid stance</li><li>• User feels they are walking up hill or forefoot feels excessively long</li></ul>	<ul style="list-style-type: none"><li>• Increase plantar flexion resistance</li><li>• Check A-P shift alignment; ensure foot is not too anteriorly positioned</li><li>• Check distribution of plantar and dorsiflexion movement; ensure that the plantar flexion range is not excessive</li><li>• Check spring category is not too soft, if so fit a higher rate spring</li></ul>
2	<ul style="list-style-type: none"><li>• Progression from heel strike to mid stance is too rapid</li><li>• Difficulty in controlling the energy return from the foot at the heel strike (reduced knee stability)</li><li>• User feels heel is too hard, fore foot is too short</li></ul>	<ul style="list-style-type: none"><li>• Reduce plantar flexion resistance</li><li>• Check A-P shift alignment; ensure foot is not too posteriorly positioned</li><li>• Check distribution of plantar and dorsiflexion movement; ensure that there is adequate plantarflexion range</li><li>• Check the spring category is not too high for the weight and activity of the patient, if so fit lower rate spring</li></ul>
3	<p>Heel contact and progression feel OK but:</p> <ul style="list-style-type: none"><li>• Forefoot feels too soft</li><li>• Forefoot feels too short</li><li>• User feels they are walking down hill, possibly with reduced knee stability</li><li>• Lack of energy return</li></ul>	<ul style="list-style-type: none"><li>• Increase dorsiflexion resistance</li><li>• Check A-P shift alignment; ensure foot is not too posteriorly positioned</li><li>• Check distribution of plantar and dorsiflexion movement; ensure that there is not excessive dorsiflexion range</li><li>• Check the spring category is not too soft for the weight and activity of the patient, if so fit higher rate spring</li></ul>

	Symptoms:	Remedy
4	<ul style="list-style-type: none"> <li>• Forefoot feels too rigid</li> <li>• Forefoot feels too long</li> <li>• Feels like walking up hill</li> </ul>	<ul style="list-style-type: none"> <li>• Reduce dorsiflexion resistance</li> <li>• Check A-P shift alignment; ensure foot is not too anteriorly positioned</li> <li>• Check distribution of plantar and dorsiflexion movement; ensure that there is sufficient dorsiflexion range</li> <li>• Check the spring category is not too rigid for the weight and activity of the patient, if so fit lower rate spring</li> </ul>

## 10 Maintenance Instructions

 **Be aware of finger trap hazard at all times.**

### Footshell removal

# 10 Maintenance Instructions (cont.)

## Spring replacement

**3** Remove heel spring and screws.

**4** Remove toe spring screw, replace toe. Upon reassembly, use Loctite 243 (926012) and torque to 35 Nm. Ensure toe spring is central to the carrier.

**5** Re-assemble with replacement heel spring. Use Loctite 243 (926012) and torque to 15 Nm.

**6** Cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

**7** If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

**8** Lubricate toe and heel if required. (Footshell is pre-lubricated).


# 10 Maintenance Instructions (cont.)

**9** Slide carrier/heel spring assembly into the foot shell.

**10** Toe spring location in foot shell

**11** Use a suitable lever to encourage the heel spring into location in the foot shell.

**12** Ensure heel spring is engaged into slot.

Heel spring location slot

**13**

# 11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight [Size 26N]:	1,2 kg ( 2 lb 10 ozs)
Recommended Activity Level:	3
Maximum User Weight:	125 kg (275 lb)*
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Hydraulic Ankle Motion (excludes additional range of motion provided by heel and toe springs)	6 degrees plantar flexion to 3 degrees dorsiflexion
Build Height:	[sizes 22-26] 170 mm
[See diagram below]	[sizes 27-30] 175 mm

**Battery:**

Battery Type Rechargeable

Battery Life Variable , but at least 24hrs from full charge. Charge daily


**Battery Charger:**

Charging time to full charge 2 hours

Input Voltage 100V to 240V AC  
 Frequency 50 to 60Hz

**Fitting Length**

\* for weights above 125 kg (275 lb)  
 up to 150 kg (330 lb) contact a  
 Blatchford representative.

Sizes	A	B
22-26	170 mm	65 mm
27-28	175 mm	70 mm
29-30		75 mm

## 12 Replacement Parts

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22–24	25–26	27–28	29–30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S


Foot Shell (for dark add 'D')			Item	Part. No.
Size/Side	Narrow	Wide		
22L	539038S	-	Glide Sock (Sizes 22–26)	531011
22R	539039S	-	Glide Sock (Sizes 27–30)	532811
23L	539040S	-		
23R	539041S	-	Bluetooth® Dongle [Micro-D]	409093
24L	539042S	-	Micro-D Connector Cover	406433
24R	539043S	-		
25L	539044SN	539044SW	Battery Charger Kit [Micro-D]:	409087E
25R	539045SN	539045SW	–Battery Charger	950215
26L	539046SN	539046SW	–Charger Cable Adapter [Micro-D]	409091
26R	539047SN	539047SW	–Battery Charger Adaptor (UK)	950103
27L	539048SN	539048SW	–Battery Charger Adaptor (USA)	950104
27R	539049SN	539049SW	–Battery Charger Adaptor (EURO)	950105
28L	-	539050S		
28R	-	539051S	Elan Comms Kit [Micro-D]:	406431E
29L	-	539052S	–Comms Cable (USB to TTL)	950928
29R	-	539053S	–Comms Cable Adapter	409078
30L	-	539054S	–USB Bluetooth® Adapter	950129
30R	-	539055S	–Interface Cable [Micro-D]	409084
			–USB Key - Loaded	406413E

# 13 Installing and using the Elan Comms Cable

## Introduction:

The Comms Cable can be used for setting up the Elan via a USB port with a laptop computer with no Bluetooth capability or where it is not possible to make a Bluetooth connection.

Before using the cable install the driver file from the Elan USB software stick.

## To install the Comms Cable software on a laptop computer

**DO NOT** use a laptop computer or similar device with an external power supply when connected to the Elan with the comms cable.

Before installing the Comms Cable Driver Software.

- DO NOT** connect the Comms Cable to the computer
- DO NOT** launch the Elan Software Interface.

### 1 To install the Comms Cable Driver file

Put the Elan USB software stick into a spare slot in the Laptop.

In the Folder *EN*

double-click on the Comms Cable driver setup file:

**CDM20814\_Setup.exe.**

### 2 Connect the laptop to the Elan using the Comms Cable(USB), Comms Cable Adapter and Micro-D Interface Cable.

### 3 In Windows Device Manager

Windows 8:

Settings > PC Info > Device Manager > Ports (COM & LPT)


**Windows 7:**

Start > Computer>System Properties> Device Manager>Ports (COM & LPT)

Identify and make a note of the

**USB Serial Port number** e.g COM47

The Comms Cable is now installed.

## 13.1 Using the Elan PC Software Interface

**DO NOT** use a laptop computer or similar device with an external power supply when connected to the Elan.

**Take all necessary precautions when moving and walking with the Elan connected by cable to a laptop computer to avoid accident or injury.**


**Note:**  
Disconnecting the comms cable from the Elan user

(Software versions 1.0.5 or later)

It is possible to disconnect/reconnect the Elan from the laptop and comms cable during set up, except when the calibration sequence is updating the settings and the update Elan settings dialogue box appears (See Section 9.5 for the calibration sequence). Disconnect when the *calibration ready* rotating symbol appears, reconnect to update and commence the next stage.

Only disconnect/reconnect from the phono jack at the Elan not from the USB. Disconnecting at the USB will lose the Com port connection requiring a new log in.


Updating calibration settings dialogue box

calibration ready rotating symbol

### Launching the Interface

- 1 Plug the Comms Cable, adapter and Micro-D Interface Cable into the Elan and a USB port on the laptop then launch the Software Interface Program on the laptop computer (see Section 9.1).
- 2 At the *Login to Elan Interface* screen Enter Name and Password if required then Click on [*Login*]
- 3 At the *Bluetooth Connection* screen Click on [*Auto connect*] **A**  
or  
Select the previously noted **USB Serial Port Number (COM\_)** from the pull down Menu. **B**  
The visual indicators at the bottom of the screen will be active on connection.
- 4 Continue the set-up and calibrate the Elan for walking see Sections 9.3-9.6.


visual indicators

# 14 Options

Click on:

to select the *Options* screen.

## 1 Battery Indicator

Shows battery charge level with a flashing blue LED when foot is plantar flexed.

## 2 Standing Support Mode

Default is OFF

Increases resistance to reduce fatigue when standing for long periods.

Standing mode is automatically disabled when battery charge becomes low to conserve remaining power.

## 3 Fast Walk Mode (TT users only)

Enhances resistance settings to assist fast walking.

## 4 Ramp Assist Mode

Enhances resistance settings to assist ascending a ramp/incline.

## 5 Ramp Brake Mode

Enhances resistance settings with a braking effect when descending a ramp/incline.

## Mode default settings

Mode	TT	TF
Battery	ON	ON
Standing Support	OFF	OFF
Fast Walk	ON	n/a
Ramp Assist	ON	ON
Ramp Brake	ON	ON

### Note:

Disabling [OFF] the Fast Walk and Ramp Modes will disable the response to certain conditions.

Fast Walk Mode is not recommended for TF users.


## 15 Walking Calibration: earlier versions

Click on:

Calibrate

to select the *Walking Calibration* screen.

**i** Calibration may be repeated at any time to fine tune the set-up and is recommended after any change in alignment and/or basic resistance settings.

### Walking Calibration procedure:

Move the slider to select the amputation level,

**B** Transfemoral *TF* or Transtibial *TT*

#### 1 Start Calibration (*stage 1*)

To calibrate the Elan to the user's individual gait characteristics instruct the user to walk on a horizontal and level surface and at a comfortable pace.

After several steps, click on:

**C** Press to start calibration

The step counter counts down the number of remaining steps required. **D**

(Note: 8-10 uninterrupted, regular steps are required for calibration).

*Continue to Stage 2*

**i** The Elan software adapts to the individual user's gait characteristics from the Calibration stages. (See Section 3 for more information)

## Walking Calibration procedure (continued)

### Walking Calibration

Walking calibration must be carried out on a horizontal and even surface at a normal walking pace

Select amputation level

TF

TT

Steps remaining: n/a

Press to start calibration

Calibration stage 2 (assist) **E**

Calibration stage 3 (brake) **F**

Factory Reset

**D** **H** 

### 2 Calibration stage 2 (assist)

Instruct the user to walk on a horizontal and level surface and at a comfortable pace and after several steps,

click on:

**E** Calibration stage 2 (assist)

The steps counter shows the number of steps remaining. **D**

### 3 Calibration stage 3 (brake)

Instruct the user to walk on a horizontal and level surface and at a comfortable pace and after several steps,

click on:

**F** Calibration stage 3 (brake)

This completes the Elan calibration.

Click on:

[Logout] and [X] to exit the Program.

 Stages 2 and 3 will feel unfamiliar; this is normal and will pass when completed.

#### Factory Reset:

Should it be necessary to begin a completely new set-up for the Elan,

#### Caution:

Clicking on [**Factory Reset**] will **delete** all the Dynamic Alignment adjustment and any walking calibration data stored in the Elan.

**H** Click on **Factory Reset**

If this option is chosen return to Section 9.4 to begin set-up.

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by the component combinations that were not authorized by the manufacturer.

## CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

## Warranty

The Elan is warranted for – 36 months – foot shell 12 months - glide sock 3 months.

See Blatchford website for the current full warranty statement.

Scheduled service every 20 months is required, for maintenance of the warranty.

For extended warranty contact service provider.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions.

This device complies with part 15 of the FCC Rules.

Operation is subject to the following three conditions:-

- 1 This device may not cause harmful interference.
- 2 This device must accept any interference received, including interference that may cause undesired operation.
- 3 Contains pre-approved Bluetooth® module FCC ID: T9J-RN41-1

## Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

Bluetooth® is the registered trademark of Bluetooth® SIG

Elan and Blatchford are registered trademarks of Blatchford Products Limited

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© Blatchford Products Limited 2019. All rights reserved.
Endolite India Ltd