KX06

Instructions for Use

KX06
KX06-Swing

EN	Instructions for Use	3
FR	Instructions D'Utilisation	12
DE	Gebrauchsanweisung	21
IT	Istruzioni per L'Uso	30
ES	Instrucciones de Uso	39

Blatchford:

1 Description and Purpose

EN

These instructions are for the practitioner.

The term device is used throughout these instructions for use to refer to KX06.

- This device is to be used exclusively as part of a lower limb prosthesis.
- This device provides swing and (except KX06-Swing) stance control for active hip disarticulation, trans-femoral and knee disarticulation users.
- The geometry of the device provides increased toe clearance during swing phase. The build height is short and upon sitting minimises protrusion of the prosthetic knee beyond the contra-lateral limb. This combination of features reduces asymmetry for users with longer residual limbs.
- This device is designed for up to Activity Level 4 use.
- This device is manufactured from water resistant materials and is designed to accommodate limited exposure to such environments.

Contra-indications

- Residual muscular weakness, contractures or proprioceptive dysfunction including poor balance.
- Contra lateral joint instabilities or pathology
- Complicated conditions involving multiple disabilities
- Activity Level 1 users

Ensure that the user has understood any instructions for use, drawing particular attention to the safety information.


Product code

- | | |
|--------------|---|
| • KX06 | Fitted with Swing and Stance Knee Control |
| • KX06-Swing | Fitted with Swing Only Control |

2 Construction

Principal Parts:

- Frame Aluminum Alloy, Brass, Stainless Steel, Titanium, Steel
- Knee control Various materials principally Aluminum Alloy
Stainless Steel, Acetal Homopolymer, NBR, Hydraulic Fluid
- Knee pad Urethane

Cosmesis

The following standard Blatchford range of continuous polyester foam cosmesis may be used:

- 561021 - Fairing Unshaped Long
- 561022 - Fairing Unshaped Short
- 561023 - Fairing Unshaped Slim Long

3 Function

- The device has a polycentric design comprising pivots connecting a chassis to a main body via 4 link arms. A hydraulic knee control sits within this frame assembly, fitted between the rear proximal axis and a pivot distally within the anterior link arms above the lower anterior pivot.
- The geometry of the device provides efficient progression into swing phase at the required point in the gait cycle, 18mm of increased toe clearance during the swing phase when compared to a single axis knee and short protrusion when sitting.
- The device has dual stance control combining the polycentric geometrical variable instantaneous centre of rotation with an adjustable hydraulic yielding support, initiated at knee extension prior to heel strike and released by an extension moment at toe off. The KX06-Swing version provides only geometric stability during stance.
- Both versions of the device provide an adjustable hydraulic swing phase control to accommodate variable cadence.

4 Safety Information

The caution symbol highlights safety information which must be followed carefully.

Be aware of finger trap hazard at all times.

Any changes in performance of the device e.g. instability or double action in the device, should be immediately reported to the practitioner.

Always use a hand rail when descending stairs and at any other time if available.

Any excessive changes in heel height may adversely affect the stability of the device.

Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.

The user should be advised to contact their practitioner if their condition changes.

5 Maintenance

Maintenance must be carried out by competent personnel.

A visual inspection annually is recommended.

Check for visual defects that may affect proper function.

A loaner system is available should servicing be required.

The user should be advised:

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in stiffness
- Instability
- Any unusual noises

Cleaning:

Use a damp cloth and mild soap to clean the outside surfaces.

DO NOT use aggressive cleansers.

If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried.

6 Limitations on Use

Intended Life:

- Service life of this device is covered by the warranty period.
- This device is recommended for use with other Blatchford products.
- A local risk assessment should be carried out based upon activity and usage.

Lifting Loads:

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment:

Avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Corrosion:

This device has been tested to ASTM B117-07a Neutral Salt Spray Test (BS EN ISO 9227:2006).

This device is not designed for prolonged submersion. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Light surface corrosion will not affect the function or security of this device.

However, if heavy corrosion is evident, stop using this device and contact your practitioner.

Suitable for shower use

Operating and Storage Temperature Range:

Exclusively for use between temperatures of -10°C and 50°C (14°F and 122°F)

7 Alignment and Setup

Users be aware of potential finger trap hazard.



Note:

Refer to fitting instructions supplied with the hydraulic cylinder for adjustment of Stance/Swing settings.

Proximal attachment

When fitting a proximal alignment component using the 4-bolt option, ensure the securing screws protrude by a maximum of 6 mm and a minimum of 5 mm to ensure correct thread engagement into the chassis. Two short screws are supplied with the device if required for the two posterior holes in the chassis; for example when using the 4-bolt pyramid adaptors P189310, 189310-C or 189310-R.


Discard screws if not required

Note:

If using an existing socket, that does not take this 22 mm dimension into account, an A-P slide adaptor may be necessary to achieve the desired TKA line.

Distal pyramid attachment

Note:

The dome is loose until the pyramid screws are tightened. Ensure anti-rotation 'peg'* is engaged in the slot upon tightening.

8 Cylinder Replacement

Note:

Ensure the posterior links are supported upon removal and replacement of the pivot pin 235420 to prevent misalignment. The pin is a slight interference fit in the bearings. Complete disassembly of this pivot is not required, remove one screw and ease pin only far enough through assembly to release cylinder before replacement.

Special tools:


- 940084 Upper Pivot Alignment Pin
- 940274 Torx T20 (3/8" Drive)
- 941252 Clamp blocks (optional)

CYLINDER

932281 Swing & Stance

or

932282 Swing only


9 Socket & Knee Pad

Socket attachment

Ensure socket does not contact hydraulic cylinder or links when fully flexed.


The socket may contact the body only.

Knee pad replacement

To remove Knee pad remove the 2 securing screws as shown.

When refitting avoid over tightening.


10 Technical Data

Operating and Storage Temperature Range:	-10°C to 50°C 14°F to 122°F
Component Weight:	1.23 kg (2lb 10oz)
Activity Level:	3–4
Maximum User Weight:	
Activity Level 3	150 kg (330 lb)
Activity Level 4	145 kg (320 lb)
Maximum flexion angle:	160 degrees
Proximal Alignment Attachment:	Chassis to accommodate single bolt and 4 bolt options
Distal Alignment Attachment:	Male Pyramid (Blatchford)
Build Height:	
Knee to Distal pyramid	230 mm

Materials:
Aluminum Alloy, Stainless Steel, Titanium, Steel, Urethane, Brass, Acetal Homopolymer, NBR and Hydraulic Fluid

Key dimensions

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CE Conformity

This product meets the requirements of 93/42/EEC guidelines for medical products. This product has been classified as a Class I Product according to the classification criteria outlined in Appendix IX of the guidelines. The Declaration of Conformity was therefore created by Blatchford Products Limited with sole responsibility according to Appendix VII of the guidelines.

Warranty

This device is warranted for 3 years.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions.

See Blatchford website for the current full warranty statement.

Manufacturer's Registered Address

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Blatchford SAS
Endolite India Ltd