

Custom fabricated high definition silicone cosmeses are now available with us. Endolite now offers custom fabricated cosmetic fingers and partial hand prosthesis.

All our cosmetic restorations are custom fabricated from a plan duplicate of the patient's stump to achieve proper proportions. Each cosmetic is sculpted to meet individual requirements and has an extra ordinarily realistic appearance and feel.

Our procedure require three dimensional impression of the entire surface of both affected and contralateral sides.

Our fabrication technique involves coloring the silicone material itself prior to sculpting, as opposed to surface painting.

This process allows the layering of different shades of silicone to provide translucent finish to the cosmesis similar to real human tissue.

