



E-Carbon Feet Range

Blatchford 

E-Carbon Feet Range

Blatchford's E-Carbon range combines the convenience of an everyday foot with a high level of energy return that also makes them ideal for recreational sports.

Each foot is based on a lightweight, responsive and durable design, featuring a split toe-spring that adapts and conforms to terrain, enhancing balance and confidence with every step. The footshell now has a sandal-toe design for greater choice of footwear.

Options include:

- A VT element that offers additional rotation and vertical shock absorption to reduce shear forces and loading rates at the socket interface. The user is able to move and adapt more freely, making it ideal for pivoting sports such as tennis and golf, or for hiking.
- A multi-axial ankle that offers shock absorption and eases residual limb pressure for a smooth, comfortable walking experience.

Elite Blade



Suitable for submersion

The waterproof Elite Blade is a lightweight, high performance flexible blade. It is designed to deflect shock from the residual limb, making it ideal for high impact use.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe	+ Top Housing
EB	25	L	N	5	S	Adapter

Foot example: Elite Blade, size 25 left, narrow, spring rating 5, sandal toe.



Elite Blade^{VT}



Suitable for outdoor use

Elite Blade^{VT} has been designed to offer all of the benefits of Elite Blade with an additional rotation and vertical shock absorption element.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe	VT Selection
EB	25	L	N	5	S	Male VTM1 or VTF1 VTM2 or VTF2 VTM3 or VTF3 VTM4 or VTF4

Foot example: Elite Blade^{VT}, size 25 left, narrow, spring rating 5, sandal toe, VT.



Elite²



Suitable for submersion

Elite 2 is lightweight, responsive foot with independent E-Carbon foot springs. Designed for high impact, this foot is ideal for a variety of recreational sports.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe
EL	25	L	N	5	S

Foot example: Elite2, size 25 left, narrow, spring rating 5, sandal toe.



Elite^{VT}



Suitable for outdoor use

Elite^{VT} combines the technology of Elite2 with an additional rotation and vertical shock absorption element.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe	Shock Spring
EVT	25	L	N	7	S	3

Foot example: Elite^{VT}, size 25 left, narrow, spring set 7, sandal toe, shock spring 3.



*Narrow (N) and Wide (W) available for sizes 25-27 only.

Javelin wearer, Bill, had always longed for a foot that offers great stability to help reduce the number of falls he had on the ice. When Bill was first introduced to Javelin he found it provided him with exceptional stability on the ice.



“It’s amazing, I’m able to pivot on my back foot a lot better than before.”

Bill, Javelin Wearer

Epirus



Epirus is a low profile, multi-axial prosthetic foot, which combines a compliant biomimetic ankle joint with independent e-carbon springs to offer a smooth, comfortable walking experience with excellent energy response.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe
EP	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Esprit



Esprit is a low profile, lightweight, waterproof foot that provides excellent energy response. Ideal for longer residual limbs or where additional components are required.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe
ESP	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Javelin



Javelin is an extremely lightweight prosthetic foot, which offers users good energy response over a range of low to moderate impact sports and activities.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe	+ Top Housing
JAV	25	L	N	3	S	Adapter

Example: foot size 25, left, narrow, spring rating 3, sandal toe, adapter.



*Narrow (N) and Wide (W) available for sizes 25-27 only.

+44 (0) 1256 316600 | customer.service@blatchford.co.uk
Blatchford Products Limited, Unit D Antura, Kingsland Business Park,
Basingstoke, RG24 8PZ, United Kingdom.

@blatchfordgrp | blatchford.co.uk

214069182 Iss2 12/19. Information correct at time of print.

Blatchford 